

O3ELITE SINGLE GAMMA CHART

O2 LPM	Ozone Generator Dial Setting							
	1	2	3	4	5	6	7	8
1	2	4	5	6	7	8	9	10
3/4	3	5	7	9	11	13	15	17
1/2	4	7	9	11	14	16	17	19
1/4	9	14	18	22	25	29	31	33
1/8	18	29	33	41	44	46	50	52
1/16	30	37	42	48	54	56	58	63
1/32	36	42	44	51	57	62	66	70

How to use this chart: Determine how many gammas of ozone you want to produce. Find that number in the chart. Set your oxygen to the LPM setting on the left and your O3Elite dial to the setting above.

Example: To produce 18 gamma, you can use 1/4 LPM of oxygen and turn your ozone dial to 3. To produce 50 gamma, you could use 1/8 LPM of oxygen and turn your ozone dial to 7. **Note:** Gamma = ug/ml

O3ELITE DUAL GAMMA CHART

O2 LPM	Dial One Only					Both Dials				
	1	2	3	4	5	5,1	5,2	5,3	5,4	5,5
1	1	3	5	7	8	11	13	15	17	18
3/4	2	4	7	9	10	14	16	18	20	21
1/2	4	8	11	13	15	19	22	26	30	33
1/4	7	13	18	23	26	32	37	42	47	53
1/8	12	24	33	40	45	51	60	68	70	75
1/16	21	40	49	55	59	73	81	85	90	95
1/32	30	52	65	71	78	88	95	101	106	110

How to use this chart: Determine how many gammas of ozone you want to produce. Find that number in the chart. Set your oxygen to the LPM setting on the left and your O3Elite dial to the setting above. **Example:** To produce 26 gamma, you can use 1/4 LPM of oxygen and turn dial 1 to 5. To produce 70 gamma, you could use 1/8 LPM of oxygen and turn dial 1 to 5 and dial 2 to 4. **Note:** Gamma = ug/ml