

Ozone Accessories

Ozone therapy can be used in a wide variety of ways. We have high quality accessories that can be used to expand your system to achieve practically any goal. Glassware, bags, adapters, tubing; if you need it, we can probably provide it. More accessories are available at www.promolife.com.



Promolife offers the best support in the industry, before and after the sale. We can help you find exactly what you need. Use our instructional videos to set everything up quickly. We are also available by phone, email or Skype for all your questions.



Most Popular

Ozone Water Bubbler

(PL-6) \$99.95

1000 ml - Used for bubbling ozone into water for drinking or sanitation purposes. Marked for easy measurement. The built in ozone destruct traps any excess ozone so it can't enter the room. Comes with an extra stopper for storage.



Large Humidifier

(PL-500HUM) 500ml \$124.95
(PL-1000HUM) 1000 ml \$139.95

A larger version of our ozone humidifier, in either 500ml or 1000ml. 500ml jar comes with a stand. The 1000ml does not need one.



Most Popular

Ozone Water Bubbler/Humidifier Combination

(PL-22) 500ml \$149.95
(PL-23) 1000ml \$164.95

Lets you switch between humidifying or ozonating water using a 3-way switch. Save money over buying a bubbler and humidifier separately.



Ozone Oil Bubbler

(PL-8) \$119.95

The 200ml oil bubbler is used to bubble ozone through oil before breathing it or for making small batches of ozonated oil. An ozone destruct is recommended to keep excess ozone from entering the room.



Ozone Trap

(PL-O3TRAP) \$85.00

Traps any excess oil from overflowing into a device. Ideal when making ozonated oils, this piece is included in the Ozonated Oil Making System.



Ozone Humidifier

(PL-9A) \$119.95

The 200ml ozone humidifier is used to humidify ozone before using it for insufflations. This makes the ozone less drying. An ozone destruct is recommended because it takes several minutes to saturate the water before use. Ozone through saturated water has the same ozone output as regular ozone.



IV Saline Drip System

(PL-OZOIV5) 500ml \$365.00
(PL-OZOIV10) 1000ml \$395.00

The Ozone IV System is used to ozonate saline solutions for IV drip modalities. Available in two sizes. Made to order. Please call or email with questions or to order.



Ozone Cups

(PL-10A) 1" \$59.95 (PL-10) 2" \$69.95
(PL-11) 3" \$79.95 (PL-12) 4" \$89.95

Cups let you focus ozone on one spot of the body. 1, 2, 3 and 4" cups available. 5" or 6" custom sizes are available upon request. Built-in destruct keeps ozone from entering the room.



Large Oil Bubbler

(PL-500OIL) 500ml \$124.95
(PL-1000OIL) 1000ml \$139.95

A larger version of our ozone oil bubbler, in either 500ml or 1000ml. The 500ml comes with a stand. The 1000ml does not need one.



Glass Ozone Syringes

(PL-50ML) 50ml \$105.00
(PL-100ML) 100ml \$125.00

All-glass syringes for exact ozone extraction. 50ml or 100ml. Also available with reusable insufflation bag and catheters.



Ozonated Oil Making System

(PL-17)

\$249.00

Make your own ozonated oils with the Promolife Ozone Oil Making System. Includes oil bubbler, stand, oil trap, destruct, tubing and accessories. The ozone trap is designed to catch any oil that might escape from the oil bubbler. The ozone destruct keeps excess ozone from entering the room. This is ideal because it can take a long period of time to fully saturate an oil.

**All of Promolife's
Ozone Glassware is
made with Class A
Laboratory Grade
Borosilicate Glass.**




Most Popular

Reusable Insufflation Bag

(PL-INSBAG)

\$79.95

Our 1000ml reusable ozone insufflation bag is made out of ozone-resistant Teflon and makes insufflations easy. Store ozone in the bag until you're ready to use it.



Kynar Limb Bag

(PL-KYN)

\$69.95

25" x 12". Ozone resistant bag to target ozone on an arm or leg; pre-installed Kynar fitting with an on and off position makes it easy to use.



3-way Stopcock

(PL-3WAYSC)

\$11.95

Three-way dial allows you to divert ozone from one application to another easily. Ozonate water in one jar, then flip the switch to ozonate oil in another or direct ozone to a destruct. There are many possibilities.



Teflon Limb Bag

(PL-TEF)

\$39.95

12"x24". The Ozone Resistant Teflon Bag for ozone therapy is 100% Teflon (PTFE). Teflon will not break down when in contact with ozone. Target ozone on an arm or leg easily.




Syringe Filling System

(PL-14)

\$125.95

The filling system allows you to easily fill multiple syringes without needing to turn off the ozone generator. Syringes attach easily via Luer lock or tubing.

Medical use only. Please contact us with questions.



Ozone Destruct System

(OZDEST)

\$49.95

A three-way switch allows you to divert ozone into a specific application or into the destruct to keep excess ozone out of the room. Ideal for humidified insufflations with standard humidifiers, which require several minutes to fully saturate the water before insufflating.



Full Body Ozone Bag

(O3BAG)

\$59.95

For administering ozone to the entire body. Simple to use, durable, three feet wide by seven feet long, with ozone resistant materials.



Earscope

(PL-16) Standard

\$49.95

(PL-16TEF) Teflon

\$49.95

Ideal for insufflating ozone into the ears. Available in standard or Teflon. Best when used with an ozone humidifier.

How do I order?

Visit us online or call us to place your order.



www.promolife.com

1-888-742-3404 | info@promolife.com

